
report_2idd Documentation

Release 2016-04

Pete Jemian

Sep 28, 2017

Contents

1	Contents	3
1.1	mda2idd_gui	3
1.2	mda2idd_report	6
1.3	mda2idd_summary	9
2	Indices and tables	11
	Python Module Index	13

Many sectors of the APS use the EPICS¹ *sscan* record² from *synApps*³ to collect their experimental data. The *saveData* process⁴ runs in an EPICS IOC to store the *sscan* record data into binary files using the *MDA*⁵ format.

This support provides a tool to read those MDA files and, for 1-D and 2-D scans, write the scan into ASCII text files according to the expected layout provided by a previous program (*yca_scanSee_report*, support based on the open-source *yorick* software with extensions to communicate with EPICS using Channel Access).




Fig. 1: Main window showing full summary of selected MDA file

¹ EPICS: <http://www.aps.anl.gov/epics/>

² EPICS *sscan* record: <http://www.aps.anl.gov/bcda/synApps/sscan/sscanRecord.html>

³ *synApps*: <http://www.aps.anl.gov/bcda/synApps/index.php>

⁴ more information on the *saveData* process: <http://www.aps.anl.gov/bcda/synApps/sscan/sscanDoc.html>

⁵ MDA format specification: http://www.aps.anl.gov/bcda/synApps/sscan/saveData_fileFormat.txt

1.1 mda2idd_gui

Source code documentation

1.1.1 mda2idd_gui Module

GUI for *mda2idd_report*

Objectives

Provide GUI tools to browse a file system and select MDA files. Process them with *mda2idd_report*.

Instructions

Browse to a directory containing MDA files. Select one. A summary will be shown. Choose: * File → Save (or ^S) to convert selected MDA file to an ASCII text file. * File → Convert entire Directory (^D) to convert all MDA files

For now, only **.mda* files may be browsed. ASCII text files will be written to directory: *../ASCII/* (relative to the MDA file directory)

Features

- presents file system list
- adds directory picker dialog and text entry box
- preview brief header or full summary of MDA file (^B)
- convert one selected MDA file to ASCII (^S)
- convert entire directory of MDA files to ASCII (^D)

Source Code Documentation

MainWindow([parent, start_fresh])Manage the application through the main window

class `mda2idd_gui.MainWindow` (*parent=None, start_fresh=False*)Bases: `wx._windows.Frame`

Manage the application through the main window

GetDefaultPreferencesFileName ()

return the name of the preferences file for this session

The preferences file, an XML file that contains recent program settings for specific program features, is saved in the HOME (or USERPROFILE on Windows) directory for the user account under the `.mda2idd_gui_rc.xml` file name. Here is an example from a Windows 7 system:

```
<?xml version="1.0" encoding="UTF-8"?>
<mda2idd_gui datetime="2013-03-06 13:09:17.593000" version="2013-02">
  <preferences_file>C:\Users\Pete\.mda2idd_gui_rc.xml</preferences_file>
  <written_by program="C:\Users\Pete\Documents\eclipse\mda2idd_
↪report\src\mda2idd_gui.py"/>
  <subversion id="$Id$"/>
  <window>
    <size h="1212" v="561"/>
    <position h="114" v="232"/>
    <sash pos="300"/>
  </window>
  <file_filter>*.mda</file_filter>
  <starting_directory>C:\Users\Pete\Documents\eclipse\mda2idd_report\data\mda
↪</starting_directory>
  <short_summary>False</short_summary>
  <mrud max_directories="9">
    <!--MRUD: Most-Recently-Used Directory-->
    <dir>C:\Users\Pete\Documents\eclipse\mda2idd_report\data\mda</dir>
    <dir>C:\Users\Pete\AppData\Local\Temp\mdalib</dir>
    <dir>C:\Users\Pete\Desktop\scanSee3.1\DATA</dir>
    <dir>C:\Users\Pete\Documents\eclipse\dc2mda\src</dir>
    <dir>C:\Users\Pete\Documents\eclipse\dc2mda\src opo</dir>
  </mrud>
</mda2idd_gui>
```

Items remembered between program sessions include:

- window size and position
- position of the sash that plits the file list from the summary output
- the list of most-recently-used directories (MRUD)
- the first directory to show (the last directory from which an MDA file was selected)

Note: If more than one copy of this program is run by the same user at the same time, the content of the preferences file will be that of the latest action that forced an update to the file content.

OnAbout (*event*)

show the “About” box

OnConvertAll (*event*)

selected the “ConvertAll” menu item

OnMenuItemExit (*event*)

User requested to quit the application

Parameters *event* – wxPython event object

OnMenuItemPrefs (*event*)

save the preferences to a file

OnMenuItemReportStyle (*event*)

OnMenuItemSave (*event*)

save the selected MDA file as ASCII

OnMrudItem (*event*)

handle MRUD menu items

OnSashMoved (*event*)

user moved the sash

OnSelectDirPicker (*event*)

user changed the text or browsed to a directory in the picker

OnSelectTreeCtrlItem (*event*)

user selected something in the directory list tree control

OnWindowGeometryChanged (*event*)

user changed the window size or position

appendSummaryText (*text*)

post new text to the summary TextCtrl, appending to any existing text

convertMdaDir (*path*)

convert all MDA files in a given directory

getPreferences (*start_fresh=False*)

set program preferences: default (*start_fresh*) and then optionally override from a file

listMdaFiles (*path*)

return a list of all MDA files in the *path* directory

messageDialog (*description, text*)

Present a dialog asking user to acknowledge something

Parameters

- **description** (*str*) – short description of message
- **text** (*str*) – message to be shown
- **yes_and_no** (*bool*) – if False (default), does not show a <No> button

readPreferences ()

read program prefs from a file

setCurrentDirectory (*directory*)
set the current directory

setStatusText (*text*)
post new text to the status bar

setSummaryText (*text*)
post new text to the summary TextCtrl, clearing any existing text

update_mrud (*newdir*)
MRUD: list of most-recently-used directories

update_mrud_menus ()
manage the MRUD menu items

writePreferences ()
save program prefs to a file

`mda2idd_gui.main()`
presents the GUI

1.2 mda2idd_report

Source code documentation

1.2.1 mda2idd_report Module

Generate ASCII text files from MDA files for APS station 2-ID-D

Objectives

- Replaces *yviewer*, `asciiRpt.py` and `mdaRpt.py`
- Creates a GUI similar to that of `asciiRpt.py` (aka *yviewer*)

Different than the output from *mdaAscii*, this module converts 1-D and 2-D scans stored in MDA files¹ into the text file format produced by `yca_scanSee_report` (a Yorick-based support).

Main Methods

- `report()`: converts MDA file to 1 or more ASCII text files, based on the rank
- `report_list()`: process a list of MDA files
- `summaryMda()`: text summary of a single MDA file (name, rank, datetime, ...)

Internal (but interesting) Methods

- `report_1d()`: report 1-D MDA scan data in the format for APS station 2-ID-D (called by `report()`)
- `report_2d()`: report 2-D MDA scan data in the format for APS station 2-ID-D (called by `report()`)
- `columnsToText()`: convert a list of column lists into rows of text

¹ MDA format specification: http://www.aps.anl.gov/bcda/synApps/sscan/saveData_fileFormat.txt

Dependencies

operating system

None. This software was developed on a Windows 7 system and tested on various Linux distributions (Ubuntu, mint, and RHEL Linux) and on MacOSX. It was also tested on solaris but the performance was too poor on that specific system to advocate its continued use.

```
>>> git clone https://github.com/BCDA-APS/mda2idd_report.git
```

Source Code Documentation

<code>summaryMda(mdaFileName)</code>	text summary of a single MDA file (name, rank, datetime, ...)
<code>report(mdaFileName[, allowException])</code>	converts MDA file to 1 or more ASCII text files, based on the rank
<code>report_1d(data)</code>	report 1-D MDA scan data in this format:
<code>report_2d(data)</code>	report 2-D MDA scan data in this format, one file for each detector:
<code>columnsToText(columns)</code>	convert a list of column lists into rows of text
<code>writeOutput(path, filename, output)</code>	write the output text buffer to the file
<code>getAsciiFileName(data[, detNum])</code>	return the proper text file name, based on the file name stored in the MDA data structure
<code>getAsciiPath(mdaFileName)</code>	given the path to the MDA file, return the related ASCII file path
<code>report_list(mdaFileList)</code>	process a list of MDA files

exception `mda2idd_report.RankException`

Bases: `exceptions.Exception`

this report can only handle ranks 1 and 2

exception `mda2idd_report.ReadMdaException`

Bases: `exceptions.Exception`

MDA files are all version 1.3 (+/- 0.01)

`mda2idd_report.columnsToText` (*columns*)

convert a list of column lists into rows of text

column widths will be chosen from the maximum character width of each column

Parameters `columns` (`[[str]]`) – list of column lists (all same length)

Returns `str` text block, with line separators

Example:

```
>>> columns = [ ['1A', '2A'], ['1B is long', '2B'], ['1C', '2C'] ]
>>> print columnsToText( columns )
1A 1B is long 1C
2A 2B          2C
```

`mda2idd_report.developer_test()`
 only for use in code development and testing

`mda2idd_report.getAsciiFileName(data, detNum=None)`
 return the proper text file name, based on the file name stored in the MDA data structure

Parameters

- **data** (*obj*) – MDA data structure returned by `mda.readMDA()`
- **detNum** (*int*) – (2-D only)

`mda2idd_report.getAsciiPath(mdaFileName)`
 given the path to the MDA file, return the related ASCII file path

Create the path to the ASCII directory if it does not exist. If we cannot create the ASCII dir path, return the MDA file path instead.

The default expectation is that the files are stored in this type of directory structure:

```
some/path/to/data/
  ./MDA/
    scan_0001.mda
  ./ASCII/
    scan_0001.1d.txt
```

`mda2idd_report.main()`
 handles command-line input

`mda2idd_report.report(mdaFileName, allowException=False)`
 converts MDA file to 1 or more ASCII text files, based on the rank

Parameters **mdaFileName** (*str*) – includes absolute or relative path to MDA file

Returns **dict** {mdaFileName: [asciiFileName]}

`mda2idd_report.report_1d(data)`
 report 1-D MDA scan data in this format:

```
1 ;
2 ; =====
3 ; Filename: /home/2-iddf/data12c3/Alix_NWU/mda/2iddf_0001.mda
4 ; 1D Scanno # = 1
5 ; title= (Scan # 1)
6 ; xtitle= PI_sample1_X(micron)
7 ; ytitle=
8 ; timeStamp= OCT 30, 2012 12:03:41
9 ; comment=
10 ;
11 ;
12 ; DIS:          P1          D1          D2          ...
13 ; Name:         2iddf:m38.VAL  S:SRcurrentAI  2idd:scaler1_cts1. ...
14 ; Desc:         PI_sample1_X  SR Current     ...
15 ; Unit:         micron       mA             ...
16 1             1197.47       122.374       169.800       ...
17 2             1198.97       122.347       169.400       ...
18 3             1200.47       122.318       171.000       ...
19 4             1201.97       122.289       173.600       ...
20 5             1203.47       122.777       169.000       ...
```

`mda2idd_report.report_2d(data)`
 report 2-D MDA scan data in this format, one file for each detector:

```

1 ; FILE: /home/2-iddf/data12c3/Alix_NWU/mda/2iddf_0012.mda
2 ; Title: Image#16 (2iddf:mca1.R9) - D01
3 ; Scan # = 12 , Detector # = 16 , col= 51 , row= 51
4 ; Yvalue: -1712.03 -1711.53 -1711.03 -1710.
5 ↪53 ...
6 ; Yindex 1 2 3
7 ↪4 ...
8 ; Xindex, Xvalue, Image(I,J) Array
9 ↪ ...
10 1 1289.55 0.00000 0.00000 0.00000 0.
11 ↪00000 ...
12 2 1290.05 0.00000 0.00000 0.00000 0.
13 ↪00000 ...
14 3 1290.55 0.00000 0.00000 0.00000 0.
15 ↪00000 ...

```

`mda2idd_report.report_list` (*mdaFileList*)
process a list of MDA files

`mda2idd_report.summaryMda` (*mdaFileName*)
text summary of a single MDA file (name, rank, datetime, ...)

Developed for the GUI to give the user a preview of the file before saving its data as ASCII to a text file.

`mda2idd_report.writeOutput` (*path, filename, output*)
write the output text buffer to the file

Parameters

- **path** (*str*) – absolute or relative path to directory where file should be written
- **filename** (*str*) – name of file to be written, existing file will be overwritten without warning
- **output** (*str*) – text buffer to write to file

1.3 mda2idd_summary

Source code documentation

1.3.1 mda2idd_summary Module

Generate ASCII text summaries of MDA files for APS station 2-ID-D

Source Code Documentation

<code>summaryMda</code> (<i>mdaFileName</i> [, <i>shortReport</i>])	text summary of a single MDA file (name, rank, datetime, ...)
<code>summary_list</code> (<i>mdaFileList</i>)	process a list of MDA files

`mda2idd_summary.main()`
handles command-line input

`mda2idd_summary.summaryMda (mdaFileName, shortReport=True)`
text summary of a single MDA file (name, rank, datetime, ...)

Developed for the GUI to give the user a preview of the file before saving its data as ASCII to a text file.

`mda2idd_summary.summary_list (mdaFileList)`
process a list of MDA files

CHAPTER 2

Indices and tables

- `genindex`
- `modindex`
- `search`

m

mda2idd_gui, 3
mda2idd_report, 6
mda2idd_summary, 9

A

appendSummaryText() (mda2idd_gui.MainWindow method), 5

C

columnsToText() (in module mda2idd_report), 7
convertMdaDir() (mda2idd_gui.MainWindow method), 5

D

developer_test() (in module mda2idd_report), 7

G

getAsciiFileName() (in module mda2idd_report), 8
getAsciiPath() (in module mda2idd_report), 8
GetDefaultPreferencesFileName()
(mda2idd_gui.MainWindow method), 4
getPreferences() (mda2idd_gui.MainWindow method), 5

L

listMdaFiles() (mda2idd_gui.MainWindow method), 5

M

main() (in module mda2idd_gui), 6
main() (in module mda2idd_report), 8
main() (in module mda2idd_summary), 10
MainWindow (class in mda2idd_gui), 4
mda2idd_gui (module), 3
mda2idd_report (module), 6
mda2idd_summary (module), 9
messageDialog() (mda2idd_gui.MainWindow method), 5

O

OnAbout() (mda2idd_gui.MainWindow method), 5
OnConvertAll() (mda2idd_gui.MainWindow method), 5
OnMenuFileItemExit() (mda2idd_gui.MainWindow method), 5
OnMenuFileItemPrefs() (mda2idd_gui.MainWindow method), 5

OnMenuFileItemReportStyle()
(mda2idd_gui.MainWindow method), 5

OnMenuFileItemSave() (mda2idd_gui.MainWindow method), 5

OnMrudItem() (mda2idd_gui.MainWindow method), 5

OnSashMoved() (mda2idd_gui.MainWindow method), 5

OnSelectDirPicker() (mda2idd_gui.MainWindow method), 5

OnSelectTreeCtrlItem() (mda2idd_gui.MainWindow method), 5

OnWindowGeometryChanged()
(mda2idd_gui.MainWindow method), 5

R

RankException, 7

ReadMdaException, 7

readPreferences() (mda2idd_gui.MainWindow method), 5

report() (in module mda2idd_report), 8

report_1d() (in module mda2idd_report), 8

report_2d() (in module mda2idd_report), 8

report_list() (in module mda2idd_report), 9

S

setCurrentDirectory() (mda2idd_gui.MainWindow method), 5

setStatusText() (mda2idd_gui.MainWindow method), 6

setSummaryText() (mda2idd_gui.MainWindow method), 6

summary_list() (in module mda2idd_summary), 10

summaryMda() (in module mda2idd_report), 9

summaryMda() (in module mda2idd_summary), 10

U

update_mrud() (mda2idd_gui.MainWindow method), 6

update_mrud_menus() (mda2idd_gui.MainWindow method), 6

W

writeOutput() (in module mda2idd_report), 9

writePreferences() (mda2idd_gui.MainWindow method),
6